HELMHOLTZ SPITZENFORSCHUNG FÜR GROSSE HERAUSFORDERUNGI

The HERA transverse polarimeter



Future Circular Collider Technical and Financial Feasibility Study 2d FCC Energy Calibration, Polarization and Mono-chromatisation workshop

The HERATPOL Stefan Schmitt, DESY

FCC EPOL WORKSHOP

19-30 September 2022 at CERN

https://indico.cern.ch/e/EPOL2022

remote participation possible

Outline

HELMHOLTZ SPITZENFORSCHUNG FÜ GROSSE HERAUSFORDER



- Introduction
 - The HERA collider
 - Polarization at HERA and HERA polarimetry
 - The HERA Transverse Polarimeter (TPOL)
 - Systematic limitations of the HERA polarimeters

Disclaimer:

this talk is on HERA polarimetry, but reflects my personal opinions only. I have been working on with the POL2000 group in the years 2000-2007, mainly on the transverse polarimeter

The HERA collider


DESY.

HEI MHOI T7 SPITZENFORSC

- Operated from 1992 to 2007
- Circumference 6.3 km
- Electrons or positrons colliding with protons
- Proton: 460-920 GeV, Leptons 27.6 GeV
- Peak luminosity ~7×10³¹ cm⁻²s⁻¹
- Lepton beam polarization above 60% achieved








September 2022

S.Schmitt, HERA transverse polarimeter

HELMHOLTZ SPITZENFORSCHUNG FÜR GROSSE HERAUSFORDERU

Beam polarization at HERA

- Proton-beam: unpolarized
- Lepton beam: unpolarized at injection energy (12 GeV)
- Lepton beam acquired transverse polarization at collision energy (27.5 GeV): Sokolov-Ternov effect
- Rise-time at HERA ~40 minutes (cf. duration of a fill: ~10 hours)
- Requirement: "flat" machine → compensating magnets for H1 & ZEUS solenoids






4

Side view

Longitudinal polarization for experiments

Siberian snake of 1st kind

- First experiment making use of HERA beam polarisation: HERMES (start in 1995)
- Spin rotators: longitudinal polarization in the HERMES straight section, transverse polarization in the arcs
- Luminosity upgrade 2000-2002
 - Install spin-rotator pairs around H1 and ZEUS
 - Remove compensating coils






The Polarimeters at HERA

- Three HERA polarimeters
 - Transverse polarimeter (TPOL) 1992-2007
 - Longitudinal polarimeter (LPOL) 1995-2007
 - LPOL Cavity polarimeter operation 2006-2007

HERA-I phase (1992-2000): no spin-rotators for H1 and ZEUS



Polarimetry requirements at HERA



- Machine setup for tuning beam energy and "harmonic bumps", to maximize polarization
 - Resonably fast feedback
 - Absolute scale uncertainty is less important (5-10%)

• Transverse polarimeter (HERA-I design)


- Experiments
 - Fast and reliable monitoring of polarization during data taking
 - Colliding bunches (H1,ZEUS) and all bunches (HERMES)
 - Absolute scale uncertainty better than 2%
- Transverse polarimeter (HERA-II design) and offline analysis
- Not covered in this talk, see backup slides
- Longitudinal polarimeter near HERMES
 - LPOL cavity polarimeter

Polarimetry at HERA



HEI MHOI T7

• Make use of backward Compton scattering off a laser beam



- Laser helicity is flipped regularly
- Polarization is proportional to differences between cross section data with opposite laser helicity


Compton scattering cross section

 $\frac{d\sigma}{d\Omega} \sim \Sigma_0 + S_3 (P_Y \Sigma_{2Y} \sin \phi + P_Z \Sigma_{2Z})$

- S_3 laser beam helicity
- P_{Y} transverse beam polarization
- P_{Z} longitudinal beam polarization

 $\boldsymbol{\Sigma}_{0},\boldsymbol{\Sigma}_{2Y},\boldsymbol{\Sigma}_{2Z}$ photon energy dependent terms


Example: $\int_{40}^{10} \int_{20}^{70} \int_{40}^{70} \int_{50}^{70} \int_{50}^$



Transverse polarimeter (TPOL) setup



- Continuous-wave laser: single photon mode (Compton scattering probability per bunch <1%)
- Vertical crossing angle 3.1mrad
- Electron and photon beams are separated by dipoles
- Photon calorimeter is 65 meter away from interaction point (lead housing)
- Laser beam-dump with optical diagnostics (measure residual linear light polarization)



Electron beam Twiss parameters at IP are chosen to give small vertical beam size of photon beam at calorimeter $\sigma_v \sim 0.5 \text{ mm}, \sigma_x \sim 2 \text{ mm}$

S.Schmitt, HERA transverse polarimeter


The photon calorimeter

- Transverse beam polarisation causes spatial asymmetry in cross section (up-donw asymmetry)
- Calorimeter is split into two optically isolated halves
- Shower-sharing between up and down depends on vertical impact point (non-linear transformation)
- Left and right channels for calibration and trigger

$$\frac{d\,\sigma}{d\,\Omega} \sim \Sigma_0 + S_3(\boldsymbol{P}_{\boldsymbol{Y}} \Sigma_{2\boldsymbol{Y}} \sin \phi)$$

 S_3 laser beam helicity

- P_{Y} transverse beam polarization
- Σ_0, Σ_{2Y} photon energy dependent terms



Transverse polarimeter online data analysis

- Polarization measurement:
 - In selected energy window: get mean of up/down asymmetry for both laser helicity states (S₃=L,R)
 - Difference of means is proportional to polarization


$$P = AP \times (\langle \eta \rangle_{S_3 = L} - \langle \eta \rangle_{S_3 = R})$$

 Analyzing power depends on beam parameters and calorimeter properties

Stat.precision: ~1% per minute (all bunches averaged)

September 2022

S.Schmitt, HERA transverse polarimeter



-0.8

-0.6 -0.4 -0.2

11

0.6 0.8

spatial coordinate

04

Operation of transverse polarimeter



- Between fills: center laser on analyzer box, measure residual linear light polarisation of L and R helicity states
- Injection and ramp: keep collimator closed, protect calorimeter
- At collision energy:
 - Adjust mirrors to maximize Compton rate (luminosity)
 - Adjust calorimeter position to have beam in its center
 - Adjust HV for calorimeter calibration
 - Measure polarisation
- Autonomous operation (autopilot)


User interface: auto-pilot, main window, details


					e	0
	Session type	 not connect 	ted On	nonitor	 control 	⊖ expert
	Global DAQ status:		UP			
	Pilot control	Pilot status	Laser	Collimators	horizontal coll.	vertical coll.
	AUTOPILOT	IDLE	OFF	CLOSED	CLOSED	DISABLED
	Lightpol control	Lightpol meas	Analyserbox	Diode	Power ^D meter	Prism
File Expert Fran run control S	OFF	IDLE	IDLE	IDLE	IDLE	ZERO
Polarizati Polarizati	Mirror center	Mirror scan	Mirror 2H	Mirror 2V	Mirror 3H	Mirror 3V
80 LPOL pol	ABOX	IDLE	CENTER	CENTER	STOPPED	STOPPED
40	Event data	Cycle	Fast DAQ	Scaler	Chopper	Pockelcell
12:00	STOPPED	STOPPED	STOPPED	RESET	OPEN Set time	OFF
	Hera leptons NOBEAM Enter new comme	15:00 16:00 DAQ status NORI D HERA leptons INJEC	20 0 10 10 10 10 10 10 10 10 10 10 10 10	20.00 21.00	short II Bean F Delay 22:00 23 Upp	th 5.0 h 30.0 min h 2.5 h 15.0 min n 1.0 h 5.0 min n 30.0 min 2.5 min variables Change Y-axis position Change Y-axis position Hide axis
P(colliding)	0.0 +/- 0.0 % Polariza	Pilot control AUTC	80	avy (1/0)		polarisation Change Y-axis
P(non-colliding) P(all)	-3.4 +/- 4.0 % Five Mi	FIIOL ACTION IDLE	60		coll. bu	TPOL ratio 5 🗠 linear lightpol diode (%) unch pol. 10r unch pol. 1m
On Off Rate	20.8 6.4 kHz	Collimators CLOS	40		coll. bu	unch pol. 30r unch pol. 5m bunch currei
Beam Y Offset	0.063+/-0.004 mm	Laser 10W	20	A Carlos and a carlos	non-co	oll. bunch po
Beam Spot	0.675+/-0.006 mm	Chopper OPEN Pockels Cell OFF	0.5 20.5 40.5 60.5	80.5 100.5 120.5 140.5	160.5 180.5 200.5	date display
Linear Light	1.6+/-2.5 % -3.4 +/-		-			
Luminosity	41.9+/-0.2 %	Calibration NOMINA		ما به : ما ب	officier	$\sim \sim $

Very high efficiency >95%

Detector upgrade in 2000-2002

- Spin rotators for HERMES, H1 and ZEUS: have to measure both polarisation of colliding bunches and all bunches
- DAQ upgrade: electronics from H1 luminosity system. Sampling ADC with two independent pipelines
- Digitisation at 40 MHz, readout by dedicated 20 MHz bus (fast DAQ branch) or by VME
- Per-pulse pedestal subtraction
- Detector upgrade: silicon-strip sensor
- Goal: in-situ calibration of energy-asymmetry response






for polarisation analysis



TPOL offline analysis



- Original analyzing power was based on simulations → polarization scale accurate to 8% NIM A329 (1993), 79
- Non-linear transformation, corrections from beam emittance, IP position, ...
- HERA-II upgrade: converter plate and silicon-strip detector \rightarrow in situ calibration
- New offline-analysis based on in-situ measurement of η-y transformation and energy response
- Offline Analysis power takes into account all known corrections, e.g.: η-y transformation, beam size and position



September 2022

S.Schmitt, HERA transverse polarimeter

Summary



- HERA: lepton beam polarisation above 60% for HERA-I (above 40% for HERA-II) achieved using the Sokolov-Ternov effect
- First polarimeter in operation: transverse polarimeter (TPOL)
- CW Argon laser with 10W gave Compton interaction rate up to 50 kHz (0.5% of the 10 MHz bunch-crossing rate)
- Simple calorimeter design with two optically separated halves and four channels \rightarrow build for reliability, not for precision
- Very robust design, fully autonomous operation
- Absolute scale precision below 2% reached only after adding converter plate and silicon detector for HERA-II operation [plus years of analysis]




Backup slides

September 2022

S.Schmitt, HERA transverse polarimeter

Losses from extra spin rotators and beam-beam effects (different polarization for colliding and non-colliding bunches)

Achieved Polarisation during HERA operation






Transverse calorimeter online calibration

Compton edge

energy)

- Average over both helicities .
- Subtract laser-off background •
- Left (L) and Right (R) channel calibration •
 - Make sure L/(L+R) is at 0.5
 - Compton edge at expected L+R (~13.8 GeV)
- Up (U) and Down (D) channel calibration •
 - Ratio (U+D)/(L+R) is analyzed as a function of x=U/(U+D)
 - Extrapolate to x=0: D calibration
 - Extrapolate to x=1: U callibration
 - Cross-check: U+D Compton edge





HEI MHO

HELMHOLTZ SPITZENFORSCHUNG FÜL GROSSE HERAUSFORDER

Operational difficulties (selection)

- Main weakness: transport system of laser light over 200 metres
- Limited diagnostics and slow mirror controllers difficulties to steer the laser into the tunnel after long shutdowns
- Longer-term: somewhat limited laser stability, rather high cost (maintenance contract with company)
- Residual linear light polarisation: difficult to adjust optics (could have profited from 2nd Pockel's cell)

- Over longer periods: 167 m damage on Mirror M4 (close to beam). Not clear whether it was from radiation or because of the (more focussed) laser beam
- Potential weakness: laser steered to electron beam, could not monitor light polarisation in that position.
- Exit window → true laser polarisation inside vaccum not known




2.9 m

M4

Longitudinal polarimeter


- HERMES physics operation: need better polarimeter with precision 1-2%
 (TPOL procision ~8% at the time)
 - (TPOL precision ~8% at the time)
- Measure longitudinal polarization between spin rotators
- Pulsed laser, multi-photon mode
 - Per shot, the total energy of ~1000 photons is measured in a crystal calorimeter
- Asymmetry between two laser helicity states → beam polarization





Longitudinal polarimeter analysis

- Energy asymmetry is fairly robust against systematic effects, analyzing power is known analytically
- Experimental difficulties
 - Pedestal from synchroton radiation
 → data with non-charged laser
 - Timing and intensity jitter
 - \rightarrow fixed energy 100mJ per shot
 - $\rightarrow\,$ correction based on laser timing
 - Calorimeter linearity
 - \rightarrow crystal calorimeter, test beam




Result: HERA-II LPOL scale uncertainty 2%


hep/ex 1201.2894

S.Schmitt, HERA transverse polarimeter

HERA LPOL Fabry-Perot polarimeter

- Added Fabry-Perot cavity in the • electron beam-line near the original LPOL IP
- Cavity is driven by 0.7W Nd:YAG laser • (1064 nm), effective power in cavity ~3000 KW. Optical table in the tunnel.
- Use sampling calorimeter from original • LPOL setup to detect photons
- Read out and histogram calorimeter ٠ data at the HERA bunch crossing rate of 10.4 MHz \rightarrow quite difficult (FPGA and CPU limitations 20 years ago)





Cavity and optical table in HERA tunnel

September 2022


S.Schmitt, HERA transverse polarimeter

CALORIMETER

~ 60m


Fabry-Perot optical setup

- Two high-reflectiveness mirrors (R>0.999) .
- Optical components mounted on optical table. Mechanically decoupled from HERA vacuum vessel using a system of bellows
- Laser is frequency-locked to a cavity • resonance using an active feedback system
- Laser helicity is selected using a rotating guarter-wave plate
- Light polarization is measured behind the . second cavity mirror



HERMES


+ 100m



2700

LPOL Cavity analysis

- Measured energy spectrum receives contributions from
 - Compton scattering (BCP)
 - Bremsstrahlung (BGP)
 - Synchrotron radiation (SRP)
- For a given events there are contributions from 1,2,3,... superimposed photons
- Analytic fit extracts relative size of these components, calorimeter properties, beam polarization, etc





HELMHOLTZ SPITZENFORSCHUNG FÜR GROSSE HERAUSFORDERUNGE


LPOL cavity results

 LPOL cavity was commissioned rather late → not used for regular operation, only in dedicated runs

However, results were very good

- Fast and accurate measurement (every 20 s for groups of bunches)
- Statistical accuracy for a single bunch: 2% per minute
- Systematic scale uncertainty 0.9%





S.Schmitt, HERA transverse polarimeter

600

 $I_{p}(\mu A)$

S.Schmitt, HERA transverse polarimeter

.

- Build-up of polarisation in a flat machine (only Dipoles)
- Rise-time depends on circumference, energy, magnetic field strength along the ring
- Exact formula can be found in NIM A329 (1993), 79

$$\tau \approx 10 \ h \frac{E^5}{E^5}$$

C : circumference [km]
S : straight sections [km]
E : energy [GeV]

 $C(C-S)^{2}$

My personal estimates ... add a big grain of salt...

	Circumference	Energy	Risetime
HERA	6.3 km	27.5 GeV	~1/2 hour
LEP	27 km	45 GeV	~10 hours
FCC	90 km	45 GeV	~400 hours
FCC	90 km	150 GeV	~1 hour

